

**Maual Test Apparatus for Copper Strip  
 Corossion Test - ASTM D130 | ASTM D4048**  
 Stainless Steel Liquid (Oil) Bath for both Pressure  
 Vessel & test tube Procedures.  
**Model: CSC-D130-L-8**



*Corossion Series*

**Maual Test Apparatus for Silver Strip Corrosion  
 Test - IP227 and Related Methods.**  
 Stainless Steel Liquid (Oil) Bath  
**Model: SLVC-IP227-L-6**


**Conform to ASTM D130 | ASTM D4048**

Features:

- Stainless steel bath
- Temperature : Min 35'C & Max 110'c
- Bath Volume : Approx : From 18 litre
- Bath Stability of  $\pm 0.5^{\circ}\text{C}$
- Digital Temperature controller with  $0.1^{\circ}\text{C}$  resolution
- A float switch for low bath fluid level safety cut-off.
- Lid cover with holes, for test tubes and pressure vessels

Accessories Available :

- Corrosion pressure vessels
- Silicon carbide paper - 240/150 grit
- Copper and Silver test strips
- Glass test tubes for D130 & IP227
- ASTM corrosion code standard
- Suitable Lid cover
- Polishing Multi Strip Holding Vice.
- Cold-finger condensor for IP 227
- Glass cradle for IP 227.



Corossion Accessory Kit available

**Country:**

**Distributor:**

