

# The F250 EcoChiller



**Save Precious Water!  
Protect the Environment!**



**NEW!**

**Julabo**  
THE TEMPERATURE CONTROL COMPANY

# F250 EcoChiller

The F250 is the ideal choice for replacing tap water in cooling applications.

Due to its high efficiency, the F250 is an economic and environmentally friendly alternative to conventional tap water cooling.

Reduction of facility water costs enables a fast ROI while saving fresh water resources

## Features & Benefits

- Environmentally friendly operation with Low energy consumption
- No side vents, instruments can be placed right next to other equipment
- All wetted parts made of stainless steel or high grade plastic
- Compact design and small footprint
- Splash-proof keypad
- Large, bright LED display
- Easy filling
- Drain tap

## Technical Specifications

JULABO Order No.	JULABO Model	Working temperature range °C	Temperature stability °C	Cooling capacity Watt				Pump capacity Flow rate / Pressure l/min bar	Filling volume liters	Dimensions W x L x H cm	Weight (net) kg	
				+20	+15	+10	+5 °C					
9 620 025	F250	+5 ... +40	±0.5	250	240	220	210	15	0.35	1.7 ... 2.6	24 x 40 x 52	27

JULABO Labortechnik GmbH  
Eisenbahnstrasse 45  
77960 Seelbach  
Germany

Phone +49 7823 51 0  
Fax +49 7823 24 91  
info@julabo.de  
www.julabo.de

**Julabo**  
THE TEMPERATURE CONTROL COMPANY