

## Scientific Glassware

### ADAPTORS

#### ADAPTORS, REDUCTION



Adaptors made from Borosilicate glass with a ground joint finish. With a Dissimilar joint size. Other combinations can be made to order.

Cat. No.	Cone Size	Socket Size
M-RA-002	19/26	14/23
M-RA-003	24/29	14/23
M-RA-004	29/32	14/23
M-RA-005	24/29	19/26
M-RA-006	29/32	19/26
M-RA-007	34/35	19/26

#### ADAPTORS, EXPANSION



Borosilicate glass expansion adaptors. Dissimilar glass ground joints. Other combinations can be made to order.

Cat. No.	Cone Size	Socket Size
M-EA-002	19/26	14/23
M-EA-003	24/29	14/23
M-EA-004	24/29	19/26
M-EA-005	29/32	14/23
M-EA-006	29/32	19/26
M-EA-007	29/32	24/29