



## Scientific Glassware

### CONDENSERS


Made from Borosilicate glass, Condensers have a great cooling area with high efficiency. Available in a great selection of shapes and sizes from air condensers to double surface area. We will ensure that you get the one most suited for you. With our in house GLASSBLOWERS we can offer custom shapes and sizes. Send us your sketch and let us do the rest.



#### LIEBIG, CONDENSER

Manufactured from borosilicate heat resistant glass these condensers are water cooled and will come supplied with detachable screw-thread connectors for safe and easy attachment of tubing. Other sizes available on request.

Cat. No.	Approx Surface Area X10-2m2	Socket Size	Cone Size	Approx Effective Length mm	Approx Overall Length mm
M-LC-004	1.30	19/26	19/26	250	380
M-LC-005	1.00	24/29	24/29	250	390
M-LC-006	1.30	29/32	29/32	250	400
M-LC-008	2.00	24/29	24/29	395	535



#### GRAHAM, CONDENSER

High efficiency borosilicate glass coil condenser supplied complete with screw-thread connectors for safe and easy attachment of tubing. Others sizes available on request.

Cat. No.	Approx Surface Area X10-2m2	Socket Size	Cone Size	Approx Effective Length mm	Approx Overall Length mm
M-GC-02	2.0	24/29	24/29	165	310
M-GC-03	2.0	29/32	29/32	165	315
M-GC-05	4.0	24/29	24/29	320	465
M-GC-06	4.0	29/32	29/32	320	470