


Scientific Glassware

JOINTED FLASKS

We are pleased to be able to manufacture to your exact specification flasks with any configuration and number of necks including a mixture of joints, flat flanges, taps and many other options. Single neck flasks are available. However, for multi-neck or special flasks simply contact us.



FLASKS, ROUND BOTTOM, SHORT NECK



Manufactured from high quality heat resistant borosilicate glass with capacities and neck sizes indicated below. Here are just some of our flasks that we offer. For other capacities and joint sizes please [contact us](#).

Cat. No.	Flask Capacity ml	Joint Size	Approx Height mm
M-RBF-009	50	19/26	90
M-RBF-014	100	24/29	110
M-RBF-018	150	24/29	115
M-RBF-023	250	24/29	130
M-RBF-027	500	24/29	150
M-RBF-034	1000	34/35	196
M-RBF-045	5000	45/40	285

FLASKS, ROUND BOTTOM, SHORT NECK



Manufactured from high quality heat resistant borosilicate glass with capacities and neck sizes indicated below. Here are just some of our flasks that we offer. For other capacities and joint sizes please [contact us](#).

Cat. No.	Flask Capacity ml	Joint Size	Approx Height mm
M-FBF-003	100	19/26	95
M-FBF-008	150	24/29	110
M-FBF-013	250	29/32	125
M-FBF-018	500	34/35	155
M-FBF-021	1000	34/35	175

Please enquire for other sizes and combinations available to order