

Scientific Glassware

JOINTED SEPARATING FUNNELS


Funnels pressure equalising, pear shape, conical, cylindrical, are all shapes that we can offer, with a large range of taps G/P roto flow to IST PTFE straight through stopcocks.

The chart below shows a selection of our most popular funnels as we have our very own "GLASSBLOWERS" we can offer custom funnels of any shape size. We can also offer graduated lines on any of our funnels. If you can draw it then we can make it.



SEPARATING, CYLINDRICAL, FUNNELS

| Cat. No. | Capacity ml | Socket | Stopcock Type | Stopcock Boar mm |
|-----------|-------------|--------|---------------|---------------------|
| M-SCF-008 | 100 | 19/26 | IST | 3 |
| M-SCF-009 | 250 | 19/26 | IST | 3 |
| M-SCF-010 | 500 | 24/29 | IST | 3 |



SEPARATING, PEAR SHAPE, FUNNELS


| Cat. No. | Capacity ml | Socket | Stopcock Type | Stopcock Boar mm |
|-----------|-------------|--------|---------------|---------------------|
| M-PSF-003 | 250 | 19/26 | GP rotoflo | 3 |
| M-PSF-005 | 500 | 24/29 | GP rotoflo | 3 |
| M-PSF-007 | 1000 | 29/32 | GP rotoflo | 3 |

JOINTED SEPARATING FUNNELS



SEPARATING, CONICAL, FUNNELS

| Cat. No. | Capacity ml | Socket | Stopcock Type | Stopcock Bore mm |
|-----------|-------------|--------|---------------|---------------------|
| M-FSC-006 | 250 | 24/29 | IST | 3 |
| M-FSC-007 | 500 | 24/29 | IST | 3 |
| M-FSC-008 | 1000 | 24/29 | IST | 3 |



PRESSURE, EQUALISING, FUNNELS

| Cat. No. | Capacity ml | Socket | Stopcock Type | Cone Size |
|-----------|-------------|--------|---------------|-----------|
| M-PEF-009 | 500 | 24/29 | IST-3 | 24/29 |
| M-PEF-011 | 1000 | 29/32 | IST-6 | 29/32 |
| M-PEF-012 | 2000 | 29/32 | IST-6 | 29/32 |

Please enquire for other sizes and combinations available to order